

ASKIR C30 CA-MI SUCTION UNIT POLY 2X2L CDE/FOOT JARS



Ref. BCASPI34



 30/07/2025 05:53:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr

ASKIR C30 CA-MI SUCTION UNIT POLY 2X2L CDE/FOOT JARS



Ref. BCASPI34

Intended for use by professionals.

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040139203	HS Code	90189084
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	9.8
Bare Unit / Weight (kg)	6.2		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			

30/07/2025 05:53:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr