

WEIGHT 400 G FOR SCHERBACK HANDLE/VAGINAL SPECULUM

Ref. IP05400



30/07/2025 01:34:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

WEIGHT 400 G FOR SCHERBACK HANDLE/VAGINAL SPECULUM



Ref. IP05400

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040121321	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	9	Packed unit / Height (cm)	4
Packed unit / Length (cm)	18	Packed unit / Weight	0.077
Bare Unit / Weight (kg)	0.077		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 01:34:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr