

GREY-TURNER FORCEPS CURVED 18 CM

Ref. IP17518



30/07/2025 07:07:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

GREY-TURNER FORCEPS CURVED 18 CM



Ref. IP17518

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040134017	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	9	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	25	Packed unit / Weight	0.045
Bare Unit / Weight (kg)	0.042		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 07:07:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr