


SCHMIEDEN SCISSORS 17.5 CM

Ref. IC24417



30/07/2025 08:25:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SCHMIEDEN SCISSORS 17.5 CM



Ref. IC24417

Schmieden Scissors 17.5 cm

Description

Schmieden Scissors 17.5 cm

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040010854	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	17.6
Packed unit / Weight	0.037	Bare Unit / Weight (kg)	0.037
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 08:25:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr