

MASTER PALM 3 OXYMETRE PORTABLE NELLCOR OXIMAX INSIDE



Ref. 223001



30/07/2025 02:34:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

MASTER PALM 3 OXYMETRE PORTABLE NELLCOR OXIMAX INSIDE



Ref. 223001

The MASTER PALM 3® pulse oximeter was developed to meet the needs of healthcare professionals. Featuring Nellcoret NellcorOximax technology, it is condensed with technology.

Description

The Master Palm 3® pulse oximeter is a suitable material for healthcare professionals compatible with NELLCOR OXIMAX technology. Its technology makes it possible to intervene quickly on all types of patients from adults to neonat, in order to accurately record their oxygen saturation and their frequency of hens.

Lightness and exceptional autonomy, you will appreciate the speed and precision of its measurements combined with its versatility.

• Adult and neonatal

technology • Exceptional autonomy • HD display with auto-rotate function that optimizes space utilization and readability by providing a vertical or horizontal layout.

• 2 display modes: Numbers and SpO2 Curve • Pulse, plethysmography curve, SpO2 and SpO2

curve display • Fully configurable, it allows you to adjust the alarm, display, sound, system

settings • Monitoring function (alarm and memories)

• Lightweight and compact 58.5 x 123 x 28 mm • Comes with rigid adult SpO2 sensor and 3 AA

batteries Accessories are available in addition according to your preferences: rechargeable battery and charging pad, pediatric or neonatal sensors, flexible sensors or clip, carrying case, etc ...

30/07/2025 02:34:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

MASTER PALM 3 OXYMETRE PORTABLE NELLCOR OXIMAX INSIDE



Ref. 223001

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446027124	HS Code	9018199000
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.696
Bare Unit / Weight (kg)	0.105		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 02:34:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr