

## CARDIOMATE 3 - ELECTROCARDIOGRAPH 3 CHANNELS



Ref. 260003



30/07/2025 17:31:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CARDIOMATE 3 - ELECTROCARDIOGRAPH 3 CHANNELS



Ref. 260003

**CARDIOMATE® ECGs are highly reliable instruments.**

### Description

ECG Cardiomete 3® and 6® are robust and reliable ECG.

Featuring a 7" color screen tilt they allow viewing 12 lead simultaneously.

Simple to use, recording and printing ECG are triggered with only one button.

The Cardiomete 3® & 6® are Cardiomete equipped with the latest technology storage and communication (SD, USB card, network) that can record internally to 10 000 ECG, send them to a computer (optional software) or directly to an external printer.

The Cardiomete 3® & Cardiomete 6® are equipped with automatic interpretation and pacemaker detection.

The internal battery life of 3 hours, their lightness and carrying handle make Cardiomete 3® & Cardiomete 6® the perfect companion in the office or traveling.

- Carrying handheld.

- Tilting screen : 7-inch color LCD display (resolution: 800x480 px)

- SD card included :4GB SD card, 10,000 reports of ECG memory in XML / PDF / BMP / JPG / DICOM

- Autonomy :The operating life is 1080 ECG.

- Automatic interpretation and intervals calculation.

- Pacemaker detection.

- Network connection, USB & printer.

- Printing and rapid analysis.

- Acquisitions and display 12 simultaneous tracks.

- Alphanumeric keypad.

### Features

|                                    |               |                         |            |
|------------------------------------|---------------|-------------------------|------------|
| CE Marking                         | CE0123        | Medical device          | Class IIa  |
| Bar code                           | 3700446025472 | HS Code                 | 9018110000 |
| Warranty (months)                  | 24            | Country                 | China      |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                         |            |
| Packed unit / Width (cm)           | nc            | Bare Unit / Width (cm)  | 24.4       |
| Packed unit / Height (cm)          | nc            | Bare Unit / Height (cm) | 6.5        |
| Packed unit / Length (cm)          | nc            | Bare Unit / Length (cm) | 31         |
| Packed unit / Weight               | 4.68          | Bare Unit / Weight (kg) | 1.975      |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                         |            |
| Minimum temperature                | -20           | Max temperature         | 50         |
| Minimum humidity                   | 25            | Maximum humidity        | 85         |
| <b>OTHER FEATURES</b>              |               |                         |            |
| Packed unit / Quantity             | 1             |                         |            |

30/07/2025 17:31:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)