


## CARDIOMATE 12 - ELECTROCARDIOGRAPH 12 CHANNELS



Ref. 260012



30/07/2025 17:34:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## CARDIOMATE 12 - ELECTROCARDIOGRAPH 12 CHANNELS



Ref. 260012

**CARDIOMATE® ECGs are highly reliable instruments.**

### Description

ECG Cardiomate 12® are robust and reliable ECG. Featuring a 10" color screen tilt they allow viewing 12 lead simultaneously. Simple to use, recording and printing ECG are triggered with only one button. The Cardiomate 12® are Cardiomate equipped with the latest technology storage and communication (SD, USB card, network) that can record internally to 10 000 ECG, send them to a computer (optional software) or directly to an external printer. The Cardiomate 12® are equipped with automatic interpretation and pacemaker detection. The internal battery life of 3 hours, their lightness and carrying handle make Cardiomate 12® the perfect companion in the office or traveling.- Carrying handheld.- Tilting screen :10-inch color LCD display (resolution: 800x480 px).- SD card included :4GB SD card, 10,000 reports of ECG memory in XML / PDF / BMP / JPG / DICOM.- Autonomy :The operating life is 1080 ECG.- Automatic interpretation and intervals calculation.- Pacemaker detection.- Network connection, USB & printer.- Printing and rapid analysis.- Acquisitions and display 12 simultaneous tracks.- Alphabetic keypad.

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446025496	HS Code	9018110000
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	40
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	9.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	42
Packed unit / Weight	8.3	Bare Unit / Weight (kg)	4.2
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	50
Minimum humidity	25	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/07/2025 17:34:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)