

MODULE TEMPERATURE COVIDIEN FILAC 3000 FOR STATION SUNTECH CT40

SUNTECH

Ref. DSU551403



CE
0413

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
3 ANS
3
3 YEARS
GUARANTEE

29/07/2025 22:16:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

MODULE TEMPERATURE COVIDIEN
FILAC 3000 FOR STATION SUNTECH
CT40

Ref. DSU551403

Description

Covidien Filiac 3000 Tempereure Module for Suntech CT40 Station

Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	3700446032432	HS Code	9018.19.5500
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	23	Packed unit / Height (cm)	14
Packed unit / Length (cm)	19	Packed unit / Weight	0.612
Bare Unit / Weight (kg)	0.34		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		


29/07/2025 22:16:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr