


BABY TO TODDLER ELECTRONIC SCALE

Ref. EHANDI3

 **Holtex**

CE



 29/07/2025 23:32:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

BABY TO TODDLER ELECTRONIC SCALE



Ref. EHANDI3

4-in-1 baby and child scales: ideal for paediatricians and general practitioners

Description

Baby and child scale that allows to weigh and measure babies with a removable tray or to weigh and measure young children from 30 to 160 cm without the tray. •

3 weighing modes: classic mode, milk absorption and mom + baby weight • Large

capacity up to 100 kg: to weigh mom + baby •

Large ergonomic weighing tray area resistant and optimized to weigh infants (600(L) x 326(W) x 91(H) mm) •

4 touch keys including manual hold function

• Hygienic, Easy to clean, stable with 4 non-slip rubber feet, easy to carry (Carrying bag sold separately) •

Conversion of units possible to kg/lb and cm/ft •

Large LCD display (79 x 28 mm) with blue backlight, easy to read and automatically rotates 180° when weight is greater than 20 kg for convenient reading

• Power: 3 x 1.5V AAA batteries (not included)

• Auto On-Off function for power saving

Features

CE Marking	CE	Medical device	No
Bar code	3661040216669	HS Code	84231090
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	34	Packed unit / Height (cm)	16
Packed unit / Length (cm)	63	Packed unit / Weight	1.5
Bare Unit / Weight (kg)	1.6		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

29/07/2025 23:32:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr