

## HOLTEX FURNITURE KIT: COUCH EPOXY+STOOL+STEP STOOL WITH 2 STEPS, BLACK



Ref. OCABREV2E\_N



 30/07/2025 03:06:32 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

# HOLTEX FURNITURE KIT: COUCH EPOXY+STOOL+STEP STOOL WITH 2 STEPS, BLACK



Ref. OCABREV2E\_N

## Full Cabinet

### Description

Set composed of:- 1 sofa

black color (ref OSA003E\_N)- 1

stool (OS0007\_N)

-

1 step 2 steps (OS0052E)

The sofa with structure in square tube in white epoxy lacquered steel

-

Adjustable headrest 75 ° with 2 self-locking mechanical racks.

-Head and seat padded with very high density foam (30 kg / m3), whose thickness is 5 cm, covered with washable plasticized fabric, fireproof (fire resistance classification: M1)

-Removable roller holder (for rolls 50 cm wide) positionable at the head or feet

- Max load: 150 kg

The stool with 5 swivel casters,

-structure 5 branches ABS, 58cm in diameter- 1

gas piston for a variable height 53 to 70 cm

-1 foam seat 30 kg / M3, diameter 43 cm, with screwed steel plate

-Plasticized coating, washable, fireproof classification M1- Max load: 150 kg

The running board 2 steps

-white epoxy lacquered steel and non-slip plastic coating

- Load: 150 kg

Stainless steel or Epoxy: What's the difference?

Epoxy: the economic choice. The steel is simply coated with a resin to delay corrosion.

Stainless steel: it is steel modified to become completely STAINLESS steel, it does not rust.

And did you know?

Today, more than 90% of stainless steel is made from recycled steel, it's good for the environment!



30/07/2025 03:06:32 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France



www.gsh-med.fr



+33(0)4 42 90 31 31



info@gsh-med.fr

## HOLTEX FURNITURE KIT: COUCH EPOXY+STOOL+STEP STOOL WITH 2 STEPS,BLACK



Ref. OCABREV2E\_N

### Features

|                             |               |                           |            |
|-----------------------------|---------------|---------------------------|------------|
| CE Marking                  | NONE          | Medical device            | No         |
| Bar code                    | 3661040216997 | HS Code                   | 9402900000 |
| Warranty (months)           | 24            | Country                   | China      |
| DIMENSIONS AND WEIGHT       |               |                           |            |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc         |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | nc         |
| Bare Unit / Weight (kg)     | 33.6          |                           |            |
| SPECIFIC STORAGE CONDITIONS |               |                           |            |
| Minimum temperature         | NA            | Max temperature           | NA         |
| Minimum humidity            | NA            | Maximum humidity          | NA         |
| OTHER FEATURES              |               |                           |            |
| Packed unit / Quantity      | 1             |                           |            |

30/07/2025 03:06:32 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)