


EXAMINATION LAMP LED 3W HOLTEX - BLACK



Ref. OSLAMP01



29/07/2025 22:15:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

EXAMINATION LAMP LED 3W HOLTEX - BLACK



Ref. OSLAMP01

Traditional and most popular medical examination lamp, with handle on the lamp holder, integrated into the switch.

Description

Light weight, easy movement, focus on any angle by flexible arm, convenient operation.

- Infrared free LED technology: it reduces the feeling of heat for the patient and the doctor. The tissues do not dry out as quickly during the procedure under very weak light radiation.

- Positioning and long range: the articulated arm ensures the ease of movement and positioning of light and thus avoids blind spots.

- Variability of lighting: the dimmer located on the control box makes it possible to reduce or accentuate the lighting power.

- Energy saving: the lifespan of an LED source is 40 times longer than that of a halogen bulb. Low power consumption: 1/3 less than conventional halogen system.

- Field of application: close examination, above the bed, intensive care, recovery rooms, dermatology, care for the elderly, gynecology, ophthalmology, ENT, minor surgery... etc.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040212814	HS Code	94054090000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	12	Packed unit / Height (cm)	111
Packed unit / Length (cm)	33	Packed unit / Weight	7
Bare Unit / Weight (kg)	3.505		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/07/2025 22:15:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr