

## EXAMINATION LAMP LED Q3 WHITE (WITH WHITE METAL STAND 5 WHEELS)



Ref. OSLAMP02



30/07/2025 10:21:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## EXAMINATION LAMP LED Q3 WHITE (WITH WHITE METAL STAND 5 WHEELS)



Ref. OSLAMP02

**Traditional and most popular medical examination lamp, with handle on the lamp holder, integrated into the switch.**

### Description

Light weight, easy movement, focus on any angle by flexible arm, convenient operation.

The articulated arm ensures the ease of movement and positioning of the light and prevents any drift. NO blind spot exams and no long reach are presented.

Our lamp provides "daylight" lighting for more accurate color rendering of fabrics and easier diagnosis and treatment.

IR-free LED technology (infrared): it reduces the feeling of heat for the patient and the doctor. The tissues do not dry out as quickly during the procedure under very weak light radiation.

- Positioning and long range: The articulated arm ensures the ease of movement and positioning of light and thus avoids blind spots.
- Variability of lighting: The dimmer located on the control box makes it possible to reduce or accentuate the lighting power.
- Energy saving: The lifespan of an LED source is 40 times longer than that of a halogen bulb. Low power consumption: 1/3 less than conventional halogen system.
- Field of application: Close examination, above the bed, intensive care, recovery rooms, dermatology, care for the elderly, gynecology, ophthalmology, ENT, minor surgery... etc.
- Light source: 3W LED- Lighting intensity: D = 30 cm 80 000 Lux, D = 50 cm 35 000

Lux

-Diameter of the light field: D 30 cm = 5 cm- Color temperature: 5000-6000K

-Color rendering index (Ra): 80%

-Temperature rise: 0 °

- Average life of LEDs: greater than 40,000 hours, longest: 50,000

-Voltage: 100V-260V 50/60Hz

-Flexible flexible arm: gooseneck length: 76 cm, no blind spot

-Fixed rod = 71 cm

-Control function: brightness control by dimmer- Diameter lamp head: 75 mm- Length lamp head:

110 mm

- Lamp head: single lamp head with on/off handle

- Power cord: length 2.5 m

- Feet with metal wheels (3.35 kg)

30/07/2025 10:21:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

# EXAMINATION LAMP LED Q3 WHITE (WITH WHITE METAL STAND 5 WHEELS)






Ref. OSLAMP02

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040212821	HS Code	9405429000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	12	Packed unit / Height (cm)	111
Packed unit / Length (cm)	33	Packed unit / Weight	7
Bare Unit / Weight (kg)	3.505		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/07/2025 10:21:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)