



AC/DC ADAPTOR FOR AIVIA



Ref. SSACCE013




 30/07/2025 08:24:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr

AC/DC ADAPTOR FOR AIVIA



Ref. SSACCE013

Description

Transformer 24 V DC 75W for cabinet AIVIA

Features

CE Marking	NONE	Medical device	No
Bar code	3661041009253	HS Code	901889010
Warranty (months)	3	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Packed unit / Height (cm)	6
Packed unit / Length (cm)	10	Packed unit / Weight	0.265
Bare Unit / Weight (kg)	0.245		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 08:24:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr