

STAINLESS STEEL MALLET DOUBLE "PLASTIC NYLON" HEAD 21 CM 200 GR



Ref. IM02400



30/07/2025 02:50:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

STAINLESS STEEL MALLET DOUBLE "PLASTIC NYLON" HEAD 21 CM 200 GR



Ref. IM02400

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040057316	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 02:50:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr