

COARCTATION CLAMP DE BAKEY 22CM ANGULATED 20°

Ref. IP19522



29/07/2025 22:22:51 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

COARCTATION CLAMP DE BAKEY 22CM ANGULATED 20°



Ref. IP19522

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040023342	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	7.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.2
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	22
Packed unit / Weight	0.064	Bare Unit / Weight (kg)	0.064
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/07/2025 22:22:51 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr