

SYRINGUES /WITHOUT ASPIRATION 2.2 CC

Ref. IDSECSL2



30/07/2025 10:19:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

SYRINGUES /WITHOUT ASPIRATION 2.2 CC



Ref. IDSECSL2

Description

Syringe cartridge side-loading

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040015965	HS Code	90189084
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	5.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	3.3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	12.4
Packed unit / Weight	0.094	Bare Unit / Weight (kg)	0.094
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 10:19:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr