

COLLIN MALLET DOUBLE HEADS 20 CM 500 GR (STAINLESS STEEL)

Ref. IM01250



30/07/2025 09:46:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

COLLIN MALLET DOUBLE HEADS 20 CM 500 GR (STAINLESS STEEL)



Ref. IM01250

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040020389	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	6	Bare Unit / Width (cm)	6.7
Packed unit / Height (cm)	4	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	21.5
Packed unit / Weight	0.637	Bare Unit / Weight (kg)	0.637
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 09:46:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr