

LITHIUM BATTERY FOR SPOT-CHECK 2VS



Ref. 546029

 30/07/2025 02:16:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LITHIUM BATTERY FOR SPOT-CHECK 2VS



Ref. 546029

Description

Replacement lithium battery for SPENGLER spot-check station 2VS

Features

CE Marking	NONE	Medical device	No
Bar code	3700446029708	HS Code	850650
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.08
Bare Unit / Weight (kg)	0.08		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 02:16:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr