

4VS MONITOR NIBP

Ref. 551600



CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
18
18 MONTHS
GUARANTEE

30/07/2025 08:39:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

4VS MONITOR NIBP



Ref. 551600

Practical for hospitals thanks to its base it is easy to move it or even position it near the patient's bed

Description

The 4VS PNI vital station has an LCD screen for better reading, visual and audible alarms are also triggered when the parameter has exceeded its limit. It offers you remarkable reliability and accuracy for nearly 10 hours, it is possible to recharge it on mains. •

- 1 Vital signs: PNI•
- 3 Blood pressure measurement modes, infants, child and adult•
- Possibility to program physiological alarms• 4 sizes of ambands: neonatal, pediatric, adult, large adult (delivered with adult cuff)
- 3500 groups of data saved•
- Battery: 10h• Size: 180 mm . 250 mm . 180 mm; Weight: 2.0 kg (with battery)
- Roller stand, wall mount, optional printer

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446032944	HS Code	9018199560
Warranty (months)	18	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	3.76
Bare Unit / Weight (kg)	2		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-40	Max temperature	55
Minimum humidity	5	Maximum humidity	95
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 08:39:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr