

LITESTICK PENLIGHT SPENGLER

Ref. 570624



 29/07/2025 22:56:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LITESTICK PENLIGHT SPENGLER



Ref. 570624

Description

This penlight is dedicated to general auscultation

- LED lighting for cold , intense and white light close to natural light
- Robust, made in aluminium
- light and compact
- on/off button
- 2 LR03 AAA batteries included

Features

CE Marking	CE	Medical device	No
Bar code	3700446009854	HS Code	901890
Warranty (months)	24		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	5	Packed unit / Height (cm)	5
Packed unit / Length (cm)	20	Packed unit / Weight	0.046
Bare Unit / Weight (kg)	0.015		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	10	Max temperature	40
OTHER FEATURES			
Colors	Gray	Packed unit / Quantity	1

29/07/2025 22:56:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr