

OXYSTART® - DIGITAL FINGERTIP PULSE OXIMETER



Ref. 221051

CE
0123

DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
2 ANS
2 YEARS
GUARANTEE



30/07/2025 06:14:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

OXYSTART® - DIGITAL FINGERTIP PULSE OXIMETER



Ref. 221051

For SpO2 spot check.

Description

For SpO2 spot checks.

Easy to use. Simply press the button to switch on.

Reading comfort measures thanks to the rotatable TFT 2 colors screen :

- Simultaneous display in less than 10 sec. :

- Pulsed oxygen saturation (% SpO2),
- Heart rate,
- The pulse strength (bar graph),
- Plethysmography curve.- Low battery indicator.

- Customizable Display: 6 display modes.

- 5 brightness levels.Switches off automatically after 5 seconds when not in use.

Comes with wrist-strap, cover and 2 batteries.

Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	3700446025465	HS Code	90181910
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	7	Packed unit / Height (cm)	4
Packed unit / Length (cm)	11	Packed unit / Weight	0.13
Bare Unit / Weight (kg)	0.06		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--

30/07/2025 06:14:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr