


## OXYSTART® - DIGITAL FINGERTIP PULSE OXIMETER BLUBERRY



Ref. 221052



30/07/2025 03:12:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## OXYSTART® - DIGITAL FINGERTIP PULSE OXIMETER BLUBERRY



Ref. 221052

For SpO2 spot check.

### Description

HIGH MEASUREMENT

ACCURACY• Plethysmography curve display•

High sensitivity sensor

EASY TO USE AND

CONVENIENT• Compact and lightweight•

Switch on at the touch of the clock•

Automatic shutdown after 5 seconds of non-use•

Automatic display

rotation• Simultaneous display in less than 10 sec. :

-pulsed oxygen saturation (% SpO2),

- heart rate,

- the intensity of the pulsatile signal (graph bar),

- the plethysmography curve.

• TFT display•

Low battery

indicator• Customizable display: 6 display modes.

• 5 levels of brightness

ADDITIONAL

INFORMATION• Comes with a wrist strap, a cover and 2 AAA LR03 batteries

### Features

|                             |               |                           |           |
|-----------------------------|---------------|---------------------------|-----------|
| CE Marking                  | CE0123        | Medical device            | Class IIb |
| Bar code                    | 3700446041809 | HS Code                   | 90181910  |
| Warranty (months)           | 24            | Country                   | China     |
| DIMENSIONS AND WEIGHT       |               |                           |           |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc        |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | 0.135     |
| Bare Unit / Weight (kg)     | 0.06          |                           |           |
| SPECIFIC STORAGE CONDITIONS |               |                           |           |
| Minimum temperature         | -40           | Max temperature           | 60        |
| Minimum humidity            | 5             | Maximum humidity          | 95        |
| OTHER FEATURES              |               |                           |           |
| Packed unit / Quantity      | 1             |                           |           |

30/07/2025 03:12:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr