

OXYFROG - DIGITAL FINGERTIP PULSE OXIMETER SPENGLER



Ref. 221056



CE
0123

DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
2 ANS
2 YEARS
GUARANTEE



30/07/2025 06:37:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

OXYFROG - DIGITAL FINGERTIP PULSE OXIMETER SPENGLER



Ref. 221056

Ideal for measuring the pulse of children

Description

High measurement accuracy:

- High sensitivity sensor
- Suitable for children from 10 kg to 40 kg

Easy to use and convenient:

- Compact & lightweight
- Automatic shut-off after 5 seconds
- Easy to use: 1 button

Additional information:

- Comes with a wrist strap and 2 AAA LR03 batteries

Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	06945040105442	HS Code	90181910
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Bare Unit / Width (cm)	3.7
Packed unit / Height (cm)	4	Bare Unit / Height (cm)	3.5
Packed unit / Length (cm)	11	Bare Unit / Length (cm)	5.9
Packed unit / Weight	0.135	Bare Unit / Weight (kg)	0.05
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 06:37:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr