

WIFI ANTENNA USB FOR ECG CARDIOMATE SPENGLER



Ref. 260013



30/07/2025 10:30:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

WIFI ANTENNA USB FOR ECG CARDIOMATE SPENGLER



Ref. 260013

=====

Description

USB WIFI antenna for ECG Cardiometre 3, 6 and 12 tracks.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446042080	HS Code	9018110000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.005		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 10:30:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr