


CARDIOMATE PC ECG SPENGLER

Ref. 260039



CE
0123



29/07/2025 22:36:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

CARDIOMATE PC ECG SPENGLER



Ref. 260039

A SPENGLER PC-ECG with XXL performance in the palm of your hand

Description

Digital ECG, the Cardiomate PC allows the real-time display of a 12-lead ECG on the screen of a computer and its printing on an A4 printer. Compatible with a Windows system, the Cardiomate PC connects via a simple USB cable. With the same performance as the Cardiomate range and with its compact and lightweight case, the Cardiomate PC is a state-of-the-art measuring instrument in the office or on the go. The software allows the visualization, recording, analysis and storage of ECGs in a clear and intuitive way.

VERY HIGH ACCURACY OF

MEASUREMENTS- Visualization of the 12 leads in real time on PC

-Sampling speed 32,000Hz

-24-bit converter

-Automatic interpretation and calculation of intervals

-Pacemaker detection

EASY TO USE

-Connection to the PC by USB

-Self-powered- Simple

and intuitive software

-Software compatible with other ECG Cardiomate

COMPACT AND LEGER

-Dimensions: 68x126x24mm

- Weight: 104g

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446042059	HS Code	9018110000
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	12.6
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2.4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	6.7
Packed unit / Weight	nc	Bare Unit / Weight (kg)	0.104
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/07/2025 22:36:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr