


## CARDIOMATE EVI 10" SPENGLER ECG TABLET (WITH PRINTER BASE)



Ref. 260051

CE  
0123



30/07/2025 03:46:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CARDIOMATE EVI 10" SPENGLER ECG TABLET (WITH PRINTER BASE)



Ref. 260051

Spengler expands its range of ECGs and offers you its latest innovation: the CARDIOMATE EVI 10" tablet, a product of a new technology, totally autonomous.

### Description

A very great simplicity of use: - 10 "color touch screen- Portable: compact size and light weight, easy to carry: 1.2 kg and 272.6 x 181.2 x 25.8 mm)- Large capacity battery, rechargeable on mains or on the recording base- Allows to print ECGs directly on a printer via wifi network or on the recorder base High-performance technology:

High accuracy: sampling speed 32,000 Hz

- Detection and marking of the

pacemaker

- Automatic interpretation (Glasgow) and calculation of intervals- Input of

patient information using a keyboard or barcode reader-

Visualization of the 12 tracks (15/18 other tablet optional) at the same time as automatic interpretation and interval calculations.

- Data export via wifi or USB port under .pdf or .xml

### Features

|                   |               |                |            |
|-------------------|---------------|----------------|------------|
| CE Marking        | CE0123        | Medical device | Class IIa  |
| Bar code          | 6971049341796 | HS Code        | 9018110000 |
| Warranty (months) | 60            | Country        | China      |

#### DIMENSIONS AND WEIGHT

|                           |     |                           |    |
|---------------------------|-----|---------------------------|----|
| Packed unit / Width (cm)  | 24  | Packed unit / Height (cm) | 36 |
| Packed unit / Length (cm) | 36  | Packed unit / Weight      | 7  |
| Bare Unit / Weight (kg)   | 3.7 |                           |    |

#### SPECIFIC STORAGE CONDITIONS

|                     |     |                  |    |
|---------------------|-----|------------------|----|
| Minimum temperature | -20 | Max temperature  | 55 |
| Minimum humidity    | 15  | Maximum humidity | 95 |

#### OTHER FEATURES

30/07/2025 03:46:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)