

DISPLAY OF 6 UNITS TENSIONIC WRIST - WRIST BPM



Ref. 527522P



30/07/2025 08:18:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

DISPLAY OF 6 UNITS TENSIONIC WRIST - WRIST BPM



Ref. 527522P

Ideal for regular self-measurement, with a very high level of precision and simple operation

Description

Display of 6 tensiometers on the wrist TENSIONIC®.
Availability 2nd quarter for the new color: Blueberry.

VERY HIGH MEASUREMENT

ACCURACY•

SPENGLER MEASUREMENT• Arrhythmia detector•

WHO tricolor risk indicator•

3 average calculations available to monitor the average of morning values

EASY TO

USE• Fully automatic•

Large backlit screen

FAMILY USE•

2 x 60 time-stamped memories

ADDITIONAL INFORMATION

• Wrist circumference 13.5-21 cm

• Comes with 2 LR03 AAA batteries and a storage case

Features

CE Marking	NONE	Medical device	No
Bar code	3700446023058	HS Code	90189010
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.815
Bare Unit / Weight (kg)	0.72		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 08:18:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr