

TENSONIC WRIST - WRIST BLOOD PRESSURE MONITOR



Ref. 527522



CE
0459



30/07/2025 05:32:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

TENSONIC WRIST - WRIST BLOOD PRESSURE MONITOR



Ref. 527522

Ideal for regular self-measurement, with a very high level of precision and simple operation

Description

For regular self-monitoring. Wrist circumference 13,5-21 cm (5,3-8,5 inch). Back-lit display. 2 x 60 memories with date and time. 3 average available : all measurements of a user / last seven daystest results for morning (5:00 am ~ 8:59 am) / last seven days test results for afternoon (18:00 pm ~ 19:59 pm). Supplied with 2 LR03 AAA batteries and a storage case.

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446022617	HS Code	90189010
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	110	Packed unit / Height (cm)	65
Packed unit / Length (cm)	165	Packed unit / Weight	0.288
Bare Unit / Weight (kg)	0.12		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	90
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 05:32:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr