


SMARTLIGHT OTOSCOPE BLACK

Ref. 570315



 30/07/2025 08:35:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SMARTLIGHT OTOSCOPE BLACK



Ref. 570315

The otoscope with conventional lighting

Description

Conventional lighting. Lightweight and robust, this otoscope is perfect for ENT auscultations especially on the move. The XHL lamp diffuses a concentrated and homogenous light beam (optional LED lamp available). Fit with an ABS and Brass head, SMARTLIGHT also comes with end clamp and built-in switch, LR6 AA batteries included.

Features

CE Marking	CE	Medical device	No
Bar code	3700446041557	HS Code	90189042
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	19	Packed unit / Height (cm)	3
Packed unit / Length (cm)	10	Packed unit / Weight	0.23
Bare Unit / Weight (kg)	0.1		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 08:35:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr