

FETAL AND VASCULAR DOPPLER WITHOUT PROBE



Ref. 743301


CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2
2 YEARS
GUARANTEE



 30/07/2025 06:15:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr

FETAL AND VASCULAR DOPPLER WITHOUT PROBE



Ref. 743301

Small and lightweight, this pocket Doppler is the essential tool for medical practices and communities

Description

For prenatal ultrasound and examination of blood vessels. Choose the probes suited to your requirements. Large colour LCD screen displaying heart rate, battery indicator, pulse wave intensity and probe frequency. Supplied with two LR6 AA batteries. Ergonomic angled probes for a comfortable use.

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446004316	HS Code	90181990
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.43
Bare Unit / Weight (kg)	0.12		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	55
Minimum humidity	5	Maximum humidity	93
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 06:15:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr