

## DIAGNOSTIC PENLIGHT LED - LUXAMED

**LUXAMED®**

Ref. D1.211.112



 30/07/2025 09:32:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## DIAGNOSTIC PENLIGHT LED - LUXAMED

**LUXAMED®**

Ref. D1.211.112

Diagnostic Penlight LED OF Luxamed is resistant and lasts over time thanks to LED technology!

This technology provides a shadow-free field of vision and optimum illumination in the area of application. The LEDs are powered by a constant current, guaranteeing considerably longer life.

### Description

Penlight for optimal eye and throat examinations, with LED light approved for eye examinations.

- Top quality ABS plastic housing
- Weight 31.5 g
- Life time LED 100.000 h
- Luminous intensity 18.000 – 22.000 mcd
- Colour temperature: approx. 7000 Kelvin
- Operating voltage 3.0 V
- Incl. 2 AAA batteries

Available in :

White (Ref : D1.211.912)

Blue (Ref : D1.211.212)

Grey (Ref : D1.211.312)

Yellow (Ref : D1.211.512)

Kids (Ref : D1.311.912)

Black (Ref : D1.211.112)

Red (Ref : D1.211.412)

Green (Ref : D1.211.612)

### Features

| CE Marking        | CE             | Medical device | Class I  |
|-------------------|----------------|----------------|----------|
| Bar code          | 04260381961673 | HS Code        | 90189084 |
| Warranty (months) | 24             | Country        | Germany  |

### Dimensions and Weight

|                          |    |                           |    |
|--------------------------|----|---------------------------|----|
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
|--------------------------|----|---------------------------|----|

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

Colors

Black