

SPIROMETER DATOSPIR AIRA D

Ref. DSSPIR01

Sibelmed



CE
0197

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
1 YEAR
GUARANTEE

29/07/2025 22:35:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SPIROMETER DATOSPIR AIRA D



Ref. DSSPIR01

For the measurement of lung flows and volumes for the diagnosis and control of respiratory diseases (asthma, COPD, etc.)

Description

- Single-use sensor (Lilly Transducer)
- Integrated W20s software with more than 40 measured parameters (CVF, CV, VMM)
- Bluetooth and USB connectivity to the computer
- Ergonomic handle for maximum comfort
- Temperature sensor integrated into the device
- Visual indicator of the status of the device
- Automatic power off, battery saving
- Comes with 25 single use transducers and 2 AAA alkaline batteries

Features

CE Marking	CE0197	Medical device	Class IIa
Bar code	3361040000038	HS Code	90181990
Warranty (months)	12	Country	Spain
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.124		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/07/2025 22:35:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr