

DISPLAY OF 12 THERMOMETRES TEMPO 10 BLUE



Ref. SP1623P



CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2
2 YEARS
GUARANTEE

30/07/2025 07:31:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

DISPLAY OF 12 THERMOMETRES TEMPO 10 BLUE



Ref. SP1623P

Description

Display of 12 TEMPO 10 blue thermometers. The TEMPO 10 thermometer

is a thermometer showing body temperature in 10 seconds•

Fever alarm•

Automatic shutdown

• Water resistant mouthpiece•

Last measurement memory•

Measurement in 10 seconds

• Comes with a storage case and button cell battery (1.5V)

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446027513	HS Code	90251920
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.65
Bare Unit / Weight (kg)	0.12		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--

30/07/2025 07:31:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr