

DISPLAY OF 12 THERMOMETER TEMPO 10 FLEX GREEN



Ref. SP1644P



CE
0123



30/07/2025 03:42:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

DISPLAY OF 12 THERMOMETER TEMPO 10 FLEX GREEN



Ref. SP1644P

Description

Display of 12 thermometers TEMPO 10 FLEX

The thermometer TEMPO 10 FLEX is a thermometer indicating the body temperature in 10 seconds and having a flexible tip:

- Fever alarm
- Automatic shutdown
- Water resistant
- Last measure memory
- 2 colors available purple or green
- Comes with: storage case, button battery (1.5V).

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446027537	HS Code	90251920
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	11	Packed unit / Height (cm)	15
Packed unit / Length (cm)	14	Packed unit / Weight	0.624
Bare Unit / Weight (kg)	0.12		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 03:42:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr