

DISPLAY OF 12 THERMOMETRES TEMPO 10 FLEX GREY



Ref. SP1648P



CE
0123



30/07/2025 01:31:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

DISPLAY OF 12 THERMOMETRES TEMPO 10 FLEX GREY



Ref. SP1648P

Description

Display of 12 thermometers TEMPO 10 Flex gray The thermometer TEMPO 10

FLEX is a thermometer indicating the body temperature in 10 seconds and having a flexible

- mouthpiece• Fever alarm•
- Automatic shutdown•
- Water-resistant and flexible tip•
- Last measurement memory
- Measurement in 10 seconds•
- 2 colors available green and gray• Comes with a storage case and button battery (1.5V)
- Comes with: Storage case, button battery (1.5V).

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446041625	HS Code	90251920
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	16
Packed unit / Length (cm)	11	Packed unit / Weight	0.66
Bare Unit / Weight (kg)	0.12		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		


 30/07/2025 01:31:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
 13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr