

DISPLAY OF 6 BLUE THERMOMETERS-TEMPO DUO

Ref. SP3100P



29/07/2025 23:42:44 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

DISPLAY OF 6 BLUE THERMOMETERS-TEMPO DUO



Ref. SP3100P

The Tempo Duo allows a frontal and ear temperature measurement by infrared emissions

Description

Display of 6 thermometers Tempo Duo

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446023133	HS Code	90251920
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.243
Bare Unit / Weight (kg)	0.3		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/07/2025 23:42:44 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr