

THERMOMETRE-TEMPO DUO BLUE

Ref. SP3101/



30/07/2025 10:22:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

THERMOMETRE-TEMPO DUO BLUE



Ref. SP3101/

The Tempo Duo allows a frontal and ear temperature measurement by infrared emissions

Description

- Infrared technology•
- Automatic calibration when switching from ear mode to temporal mode• Ear or temporal temperature reading•
- 10 memories•
- Indicator light•
- Backlit display in 3-color temporal mode according to the values: green, orange, red
- Comes with two LR03 batteries (AAA)

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446022778	HS Code	90251920
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.185
Bare Unit / Weight (kg)	0.05		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 10:22:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr