

## TEMP DUO II SPENGLER FOREHEAD AND AURICULAR THERMOMETER



Ref. SP3105



30/07/2025 07:19:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## TEMP DUO II SPENGLER FOREHEAD AND AURICULAR THERMOMETER



Ref. SP3105

Temporal or ear, the measurements are instantaneous and visible on a backlit screen.

### Description

• Measuring range: 32°C - 43°C

Equipped with a backlit screen, TEMPO DUO II turns red when the temperature is above normal. SPENGLER performance at home.

2 buttons allow you to take a quick and reliable measurement, reducing mishandling to zero. Change modes without any worries, calibration is done automatically.

With its small size and low weight, TEMPO DUO II slips into a bag for a trip. His 20 memoirs make it possible to follow his results over time. • Screen resolution: 0.1°C • Unit of measurement: C° or

F° • Memory: 20 measurements • Dimensions: 150 x

50 x 38.5 mm •

Weight: 75g with batteries •

Normal conditions of use: •

Ambient temperature: 16°C – 35°C

• Relative humidity: 85% or less

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446041786	HS Code	90251920
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	8	Packed unit / Height (cm)	5
Packed unit / Length (cm)	18	Packed unit / Weight	0.2
Bare Unit / Weight (kg)	0.06		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 07:19:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)