

DISPLAY OF 6 NON-CONTACT THERMOMETER TEMPO PRO SPENGLER



Ref. SP3110P



30/07/2025 02:46:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

DISPLAY OF 6 NON-CONTACT THERMOMETER TEMPO PRO SPENGLER



Ref. SP3110P

Description

Cardboard counter display for 6 professional thermometers TEMPO PRO Spengler

Infrared technology allows the temperature to be taken on the forehead/temple without physical contact and effectively fights against cross-contamination.

With its "auto off" mode and "low battery indicator", no unpleasant surprises during your consultations.

TECHNICAL DATA• Measurement

mode: Non-contact infrared• Measuring distance: - 5 cm• Unit of measurement: C° / F°• Display range: 32.0°C – 43.0°C• Resolution: 0.1°C• Accuracy:

0.2°C•

Memory: 30 memories•

Measurement time: Less than 1 second• Screen:

LCD

OTHER

INFORMATION• Weight: 130 g•

Dimensions: 113(W) x 51 (W) x 162 (H)•

Automatic power saving power off

• Batteries supplied with (2 x AAA batteries)

Features

CE Marking	NONE	Medical device	No
Bar code	3700446041939	HS Code	90251920
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.78		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

30/07/2025 02:46:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr