

AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR DISPLAY OF 6 UNITS

Ref. SPG340P-



30/07/2025 05:38:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR DISPLAY OF 6 UNITS



Ref. SPG340P-

SPENGLER precision in electronic version

Description

Display of 6 AUTOTENSIO® wrist model. For occasional self-monitoring. Very easy-to-use.- Quick and quiet inflation method- 120 memories with date and time- LCD screen- Fully automatic- Purple cuff- Wrist circumference 13,5 - 21 cm (5,3-8,5 inch). Supplied with 2 LR03 AAA batteries and a storage case.

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446023027	HS Code	90189010
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.7
Bare Unit / Weight (kg)	1.368		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 05:38:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr