

AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR DISPLAY OF 6 CORAL

Ref. SPG345P



30/07/2025 00:17:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR DISPLAY OF 6 CORAL



Ref. SPG345P

SPENGLER precision in electronic version

Description

Display of 6 Autotensio coral wrist.

MEASUREMENT
ACCURACY• SPENGLER MEASUREMENT•
WHO risk indicator

EASY &
CONVENIENT• Fully automatic• 2

buttons• LCD display
• 120 time-stamped memories•
Fast and quiet

ADDITIONAL INFORMATION•
Wrist circumference 13.5-21 cm
• Comes with 2 LR03 AAA batteries and a storage case

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446041601	HS Code	90189010
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.62
Bare Unit / Weight (kg)	1.368		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 00:17:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr