

AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR SPENGLER CORAL



Ref. SPG345

CE
0459

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2 ANS
2 YEARS
GUARANTEE



29/07/2025 22:46:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR SPENGLER CORAL



Ref. SPG345

SPENGLER precision in electronic version

Description

MEASUREMENT
ACCURACY• SPENGLER MEASUREMENT•
WHO risk indicator

EASY &
CONVENIENT• Fully automatic• 2

buttons• LCD display
• 120 time-stamped memories•
Fast and quiet

ADDITIONAL INFORMATION•
Wrist circumference 13.5-21 cm
• Comes with 2 LR03 AAA batteries and a storage case

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446041595	HS Code	90189010
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	10	Packed unit / Height (cm)	10
Packed unit / Length (cm)	9	Packed unit / Weight	0.27
Bare Unit / Weight (kg)	0.228		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-25	Max temperature	70
Minimum humidity	5	Maximum humidity	93
OTHER FEATURES			
Packed unit / Quantity	1		

29/07/2025 22:46:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr