

AUTOTENSIO BPM ARM CUFF ADULT (M/L) BLUE SPENGLER



Ref. SPG440-



30/07/2025 05:16:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

AUTOTENSIO BPM ARM CUFF ADULT (M/L) BLUE SPENGLER



Ref. SPG440-

Ease of use in everyday life for self-measurement family use

Description

AUTOTENSIO® upper arm model. For occasional self-monitoring. - LCD screen, - W.H.O risk indicator- Irregular heart beat detection, - Low batteries indicator, - 2 x 60 memories with date and time, - Purple cuff- Upper arm circumference 22-36 cm size M. Supplied with 4 LR06 AA batteries and a storage bag. Optional power adapter.

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446022501	HS Code	90189010
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	10
Packed unit / Length (cm)	17	Packed unit / Weight	0.594
Bare Unit / Weight (kg)	0.358		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	90
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 05:16:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr