

DISPLAY 6 UNITS AUTOTENSIO ARM-ARM BPM ADULT MEDIUM (22X36CM)

Ref. **SPG440P**



30/07/2025 04:21:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

DISPLAY 6 UNITS AUTOTENSIO ARM-ARM BPM ADULT MEDIUM (22X36CM)



Ref. SPG440P

Ease of use in everyday life for self-measurement family use

Description

Display of 6 blood pressure monitors AUTOTENSIO Blueberry.

HIGH MEASUREMENT

ACCURACY•

SPENGLER MEASUREMENT• Arrhythmia detector•

WHO tricolor risk indicator•

3 average calculations available to monitor the average of morning values EASY TO USE•

Fully automatic• 2

buttons•

New cuff 22-42 cm size M/L

FAMILY

USE• 2 x 60 time-stamped memories

ADDITIONAL INFORMATION

• Comes with 4 LR6 AA batteries and a carrying

case• Optional AC adapter

• 2 cuff sizes (S, M/L) available as an option

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	370044602034	HS Code	90189010
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	3.73
Bare Unit / Weight (kg)	2.148		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 04:21:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr