

AUTOTENSIO UPPER ARM (M/L) BLOOD PRESSURE MONITOR SPENGLER CORAL

Ref. SPG445



30/07/2025 05:55:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

AUTOTENSIO UPPER ARM (M/L) BLOOD PRESSURE MONITOR SPENGLER CORAL



Ref. SPG445

Ease of use in everyday life for self-measurement family use

Description

HIGH MEASUREMENT
ACCURACY•
SPENGLER MEASUREMENT• Arrhythmia detector•
WHO tricolor risk indicator•
3 average calculations available to monitor the average of morning values EASY TO USE•

Fully automatic• 2
buttons•
New cuff 22-42 cm size M/L

FAMILY
USE• 2 x 60 time-stamped memories

ADDITIONAL INFORMATION
• Comes with 4 LR6 AA batteries and a carrying
case• Optional AC adapter
• 2 cuff sizes (S, M/L) available as an option

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446041571	HS Code	90189010
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	10
Packed unit / Length (cm)	17	Packed unit / Weight	0.594
Bare Unit / Weight (kg)	0.358		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-25	Max temperature	70
Minimum humidity	5	Maximum humidity	93
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 05:55:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr