


Ref. C1.426.114




Rechargeable USB



29/07/2025 23:37:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr



Ref. C1.426.114

Innovative dermatoscope in pocket size, rechargeable USB

WITHOUT scale.

Description

- Illumination: LED-Ring with 6 LEDs in concentric structure for a highly detailed image, homogenous light
 - Maintenance-free, because no lamp replacement necessary
 - Removable contact plate for easy cleaning
 - Stepless zoom area with 10 x magnification for a distortion-free image
 - Optimized design of head
 - Can be dimmed
 - Automatic switchoff after 3 minutes
 - Life expectancy of LEDs: approx. 100.000 hours
- by constant current drive
- Intensity of illumination: 10.000 Lux
 - Colour temperature: approx. 4.200 K
 - Two-component handle (aluminium / fibre-glass reinforced), also combinable with LuxaScope LED Otoscope- and Ophthalmoscope head
 - Electronic ON/OFF push-button-switch
 - Power supply: 1 lithium-ion rechargeable battery 3.7 V, inclusive
 - Incl. recharger with USB plug (adapter EU/UK/US)
 - Available in white (C1.426.914) and black (C1.426.114)

Features

CE Marking	NONE	Medical device	No
Car code	04260381963431	HS Code	90189084

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	Black		