

LUXASCOPE OPHTHALMOSCOPE LED 2.5 V - LUXAMED

LUXAMED®

Ref. F1.416.954



 30/07/2025 03:47:23 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr

LUXASCOPE OPHTHALMOSCOPE LED 2.5 V - LUXAMED

LUXAMED®

Ref. F1.416.954

Direct ophthalmoscope in pocket size

Description

Direct ophthalmoscope in pocket size

Refraction balance: + 20 bis - 20 D ($\pm$ 1 2 3 4 6 8 10 15 20 D)

5 different apertures: reticle, large circle, small circle, semicircle, green filter

Illumination: LED

Can be dimmed for an individual glare-free brightness

Automatic switchoff after 3 minutes

Life expectancy of LEDs: approx. 50.000 hours by constant current drive

Colour temperature: approx. 4.200 K

Optimized design of head

Two-component handle (aluminium/fibre-glass reinforced), also combinable with LuxaScope LED otoscope- and dermatoscope head

Electronic ON/OFF push-button switch

Power supply: 2x AAA Alkaline batteries, inclusive

Available in :

- White : C1.516.914
- Blue : C1.516.934
- Pink : C1.516.974
- Yellow : C1.516.944
- Kids : C1.916.914
- Black : C1.516.114
- Green : C1.516.954
- Purple : C1.516.964

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381964407	HS Code	90185090
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Net Unit / Weight (kg)	0.12	For information purposes only, can be subject to change	

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr