

## DECONTAMINABLE ARM CUFF MONOBLOC STRAIGHT NEW BORN TOURNIQUET DDM

Ref. GBDM102



29/07/2025 23:18:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## DECONTAMINABLE ARM CUFF MONOBLOC STRAIGHT NEW BORN TOURNIQUET DDM

Ref. **GBDM102**

Monobloc straight arm cuff with one bladder and a ribbon

### Description

Monobloc PU decontaminable arm cuff are available in straight design. Made of a single monobloc unit with one bladder and a ribbon to secure the Velcro attachment.

They are available in a range of different sizes easily indentifiable by color coding of the tubing. (A : bladder width, B: bladder length, C: total length)

|                        |                | A   | B  | C  |
|------------------------|----------------|-----|----|----|
| • New born arm cuff    | White tube     | 3.5 | 29 | 38 |
| • Child arm cuff       | Pink tube      | 4.5 | 35 | 48 |
| • Small adult arm cuff | Parme tube     | 6   | 46 | 61 |
| • Adult arm cuff       | Sky-blue tube  | 8   | 46 | 67 |
| • Large adult arm cuff | Turquoise tube | 10  | 47 | 70 |

Monobloc arm cuffs can not be sterilized. They can be dry cleaned or immersed paying attention to maintain the tube connector out of water.

 29/07/2025 23:18:54 - For information purposes only, can be subject to change

## DECONTAMINABLE ARM CUFF MONOBLOC STRAIGHT NEW BORN TOURNIQUET DDM

Ref. **GBDM102**

### Features

|                             |               |                           |          |
|-----------------------------|---------------|---------------------------|----------|
| CE Marking                  | CE            | Medical device            | Class I  |
| Bar code                    | 3700468207917 | HS Code                   | 90189084 |
| Warranty (months)           | 24            | Country                   | France   |
| DIMENSIONS AND WEIGHT       |               |                           |          |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | nc       |
| SPECIFIC STORAGE CONDITIONS |               |                           |          |
| Minimum temperature         | -10           | Max temperature           | 40       |
| Minimum humidity            | 30            | Maximum humidity          | 40       |
| OTHER FEATURES              |               |                           |          |

 29/07/2025 23:18:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)