

TEMPO EAR PRO - EAR THERMOMETER

Ref. THE004



30/07/2025 03:23:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

TEMPO EAR PRO - EAR THERMOMETER



Ref. THE004

Performance, hygiene and comfort for healthcare professionals

Description

Spengler's Tempo Ear Pro is the infrared ear thermometer designed for demanding medical environments. Fast, reliable and hygienic, it makes daily temperature taking easier, whether in the office, hospital or clinic.

- Ultra-fast measurement: less than 3 seconds
- Probe preheated to 36°C: patient comfort
- Infrared technology: reliable, even in the event of movement
- Probe cover ejection: maximum hygiene without contact
- Memory of last 30 time-stamped measurements
- Automatic sleep mode after 60 seconds
- On batteries

Included with Tempo Ear Pro :

- 1 Tempo Ear Pro ear thermometer
- 1 storage base
- 20 probe covers
- 2 AA batteries
- 1 instruction manual
- Thermometer calibration and verification certificate with batch number

Also available in rechargeable "+" version, with pulse counter and safety cord

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	6970392219424	HS Code	9025199010
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	5	Packed unit / Height (cm)	3
Packed unit / Length (cm)	14	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	50
Minimum humidity	15	Maximum humidity	93

OTHER FEATURES

30/07/2025 03:23:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr