

OLSEN-HEGAR NEEDLE HOLDER 14 CM

Ref. IP32414



30/07/2025 03:48:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

OLSEN-HEGAR NEEDLE HOLDER 14 CM



Ref. IP32414

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040024936	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	14
Packed unit / Weight	0.034	Bare Unit / Weight (kg)	0.034
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 03:48:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr